

- [Investor Presentation](#)
- [Events](#)
- [Earnings Calls & Releases](#)
- [SEC Filings](#)
- [Annual Meeting & Reports](#)
- [Corporate Governance](#)
- [Investor Resources](#)


INRatio®
[Professional Products](#) [Brands](#)

- [Companies](#)
- [Professional Products](#)
- [Disease Areas](#)
- [Brands](#)
- [Acceava®](#)
- [BinaxNOW®](#)
- [Cholestech LDX®](#)
- [Clearview®](#)
- [Determine®](#)
- [DoubleCheckGold Ultra](#)
- [Immunocomb®](#)
- [INRatio®](#)
- [Inverness Medical](#)
- [TestPack®](#)
- [Matritech®](#)
- [Osteomark®](#)
- [Panbio®](#)
- [SureStep](#)
- [Triage®](#)
- [Wampole](#)
- [Distributed Products](#)
- [Health Management](#)
- [Consumer](#)

INRatio®

The top-selling INRatio® monitor is a diagnostic point-of-care system that provides prothrombin time (PT) and international normalized ratio (INR) results using fresh capillary whole blood from a fingerstick.

The INRatio® system helps oral anticoagulation management (by helping maintain the effect of anticoagulation therapy within the therapeutic range), is easy to use, and features on-board quality controls with every test. The monitor performs the PT test and two quality control tests (normal and therapeutic) simultaneously, and determines whether the controls are within preset limits. Test results are displayed in less than two minutes (as an INR only, as PT seconds and INR, or as INR, PT and QC- whichever the user prefers) thereby allowing dose adjustments face-to-face.

Recently introduced is the INRatio®2 monitor, designed for PT/INR patient self-testing at home. The INRatio®2 system allows frequent and simple tests for oral anticoagulation management, which combine ease of use with reliability of results. This means the patient and the doctor have more data and a clearer picture of PT/INR levels over time. This additional information can be used by the physician for more timely adjustments to anticoagulation therapy.

[INRatio® Monitor](#)
[INRatio®2 Monitor](#)