

- [Investor Presentation](#)
- [Events](#)
- [Earnings Calls & Releases](#)
- [SEC Filings](#)
- [Annual Meeting & Reports](#)
- [Corporate Governance](#)
- [Investor Resources](#)


ImmunoComb®
[Professional Products](#) [Brands](#)

- [Companies](#)
- [Professional Products](#)
- [Disease Areas](#)
- [Brands](#)
- [Acceava®](#)
- [BinaxNOW®](#)
- [Cholestech LDX®](#)
- [Clearview®](#)
- [Determine®](#)
- [DoubleCheckGold](#)
- [DoubleCheckGold Ultra](#)
- [ImmunoComb®](#)
- [INRatio®](#)
- [Inverness Medical](#)
- [TestPack®](#)
- [Matritech®](#)
- [Osteomark®](#)
- [Panbio®](#)
- [SureStep](#)
- [Triage®](#)
- [Wampole](#)
- [Distributed Products](#)
- [Health Management](#)
- [Consumer](#)

ImmunoComb®

Inverness Medical provides a wide range of rapid EIA (Enzyme Immunoassay) ImmunoComb® kits for the detection of antibodies or antigens in human serum or plasma.

The ImmunoComb® is a solid-phase EIA. The solid phase is a plastic comb with 12 teeth, sensitized at different spots with reactive materials and an internal control. All the reagents for the test are sealed with foil in a developing plate, containing 6 rows of 12 premixed, prefilled wells. The specimens are deposited in the wells of the first row of the developing plate.

The test is performed in a quick and easy procedure by advancing the comb from row to row at specific intervals. Results are visible within minutes, appearing as gray-blue spots on the surface of the comb's teeth. Reading the test results may be done visually or with an optical Scanner-[Combscan™](#)

- Exceptionally simple-Easy to use, no need for skilled technicians, instrument free procedure, ready to use reagents.
- Uncompromising quality- High sensitivity and specificity. Internal Control for each sample. Positive and negative Controls in each kit.
- Economic and efficient-Enables testing of small to medium batch sizes, stable results enable documentation.
- No instrumentation needed-Eliminates costly equipment.
- Flexible-Individual or multiple sample testing.
- Safe-No exposure to infectious materials.
- Efficient-Minimal contained waste.
- Broad product range-Covers the most important markers for HIV, Hepatitis, Chlamydia, ToRCH and others.

Now every lab and point-of-care unit can benefit from multiple advantages of a rapid EIA kit without the need for costly staff training procedures, and equipment.

[ImmunoComb® products](#)

[ImmunoComb® Multi-Packs](#)