

- [Investor Presentation](#)
- [Events](#)
- [Earnings Calls & Releases](#)
- [SEC Filings](#)
- [Annual Meeting & Reports](#)
- [Corporate Governance](#)
- [Investor Resources](#)


DoubleCheckGold Ultra
[Professional Products](#) [Brands](#)

- [Companies](#)
- [Professional Products](#)
- [Disease Areas](#)
- [Brands](#)
- [Acceava®](#)
- [BinaxNOW®](#)
- [Cholestech LDX®](#)
- [Clearview®](#)
- [Determine®](#)
- [DoubleCheckGold](#)
- [DoubleCheckGold Ultra](#)
- [Immunocomb®](#)
- [INRatio®](#)
- [Inverness Medical](#)
- [TestPack®](#)
- [Matritech®](#)
- [Osteomark®](#)
- [Panbio®](#)
- [SureStep](#)
- [Triage®](#)
- [Wampole](#)
- [Distributed Products](#)
- [Health Management](#)
- [Consumer](#)

DoubleCheckGold™ Ultra

DoubleCheckGold Ultra are recently introduced rapid and easy to use immunochromatography tests for simultaneous detection of HIV-1 and HIV-2 in human serum, plasma and whole blood, without instrumentation and without the need for an applicator. Highly accurate results are introduced just in 15 minutes by applying only 25 µL of sample.

The DoubleCheckGold™ Ultra HIV 1&2 offers a new level of excellence and advantages:

- Easy and rapid lateral flow assay
- Suitable for serum, plasma and whole blood samples
- Excellent performance: sensitivity 100% and specificity 99.8%
- Detection of HIV-1, HIV-2, all HIV-1 groups, subtypes (clads) and Circulating Recombinant Forms (CRF)
- Storage and test procedure at room temperature
- Results available in 15 minutes
- Simple to use and easy to handle
- Packaging of 20 tests or 100 tests

[DoubleCheckGold™ Ultra HIV 1&2](#)